



TRANSFER CHAIR **I-MOVE EZ-GO BARIATRIC**

Sturdy and Easy to handle

Adaptable Many options available

Attractive Will bring a touch of modernity and color in the unit

Motorized Very easy to use, progressive speed




I-MOVE EZ-GO is a transport chair for your patients' reception and sitting transport. It is equipped with an autonomous motorized wheel, two batteries and an external charger.

Thanks to I-MOVE EZ-GO, patients transfer is child's play ! It is particularly suitable for the transfer of overweight patients.

I-MOVE EZ-GO BARIATRIC DIMENSIONS

	LENGTH	WIDTH	HEIGHT	SEAT
I-MOVE	105 cm	90 cm	102 cm	H. 55 cm D. 46 cm W. 65 cm

PRODUCT WEIGHT : 65 Kg
SAFE WORKING LOAD: 250 Kg



Motorized wheel including battery

Our TRANSFER CHAIR range allows the reception of the patient and his transport in sitting position

Motorized transfer chair I-MOVE EZ-GO for patient welcoming

- Motorized drive system with autonomous wheel
- Switchable battery - 5 km autonomy per battery - Charge time : 5 hours
- Lithium battery without memory effect.
- Forward/Backward twist grip with progressive speed acceleration
- Steel frame epoxy coated RAL 9010
- Patient file holder
- Back wheels Ø 300 mm with central braking system, can be operated in stand position by the accompanying person
- Front swivel wheels Ø 125 mm
- Sliding and self-blocking foot rest
- Polyurethan armrests, can be folded up
- Accessory basket
- Supplied with 2 batteries and external charger
- Anti-shock protection
- Cushion in foam HR 35kg/m³ - Thickness 8 cm
- Seamless finish, vinyl coating
- Seamless finish, vinyl coating



Accessory basket



Polyurethan armrests, can be folded up



Back wheels Ø 300 mm with central braking system



Front swivel wheels Ø 125 mm
Sliding footrest

Different colors



ORANGE



GREEN




LAVANDER

Composition : Fireproof fabric jersey cotton (15%) coated with PVC (80%)
Thickness: 1,4 mm +/-0,1mm
Light resistance: >= 6



ACCESSORIES



Removable IV pole,
2-stage, 2 hook
TSR2




U-shape anti-theft, fitted with 2 hooks
infusion holder
TROAA



Holder for oxygen cylinder
TROSO2




Hooks to hold bags
PSU




Adjustable leg rest
(1 piece)
TRORJ



Wall fixation
TROMOP



Batterie supplémentaire
EZOBA



External wall charger
EZOCH



Charging station 2/3/4 chargers
(one multi-socket included)
EZOS2 - EZOS3 - EZOS4



Charging station 1 charger
EZOS1



External charger mobile version
EZOPO

REFERENCE	DESIGNATION
TSR2	Removable IV pole, 2-stage, 2 hooks
TSR4	Removable IV pole, 2-stage, 4 hooks
TROAA	U-shape anti-theft, fitted with 2 hooks infusion holder (Incompatible with TSR2)
TROSO2	Holder for oxygen cylinder
PSU	Hooks to hold bags
CLORA	Anti-static system
TROPE	Equipment for leg rest (pair)
TRORJ	Adjustable leg rest (1 piece)
TROMO - EP	Coin acceptor for 1€ coin or equivalent
TROMO - EJ	Coin acceptor for specific token
JETON	Additional specific token
TROMOP	Wall fixation
EZSBA	Without battery (depreciation)
EZOBA	Additional battery
EZOCH	External wall charger
EZOPO	External charger mobile version
EZOS1	Charging station 1 charger
EZOS2	Charging station 2 chargers (one multi-socket included)
EZOS3	Charging station with 3 chargers (one multi-socket included)
EZOS4	Charging station with 4 chargers (one multi-socket included)

A medical professional must always rely on his or her own professional clinical judgment when deciding whether to use a product(s) when treating a particular patient.

A healthcare professional anyway, refers to the package insert, product label and/or instructions for use before using any BMB MEDICAL product.

The products depicted are CE marked in accordance with EU Regulations and Directives.

The products are Class 1 Medical Devices.